



Bright, easy-to-read digital display
Masimo SET® technology
Durable metal connector

SmartMonitor® 2 Professional Series Multi-Parameter Monitoring Made Easy

SmartMonitor 2 Professional Series is the most powerful member in the Children's Medical Ventures' family of patient monitors. It can be used in the hospital or the home and is cleared for infant, pediatric and adult monitoring.

SmartMonitor® 2 Professional Series


Multi-Parameter Monitoring Made Easy

The SmartMonitor 2 PS is our top-of-the-line monitor and provides health care professionals with detailed information about a patient's heart rate, respiration and oxygen saturation. Features include Masimo SET® technology, large memory capacity, 15 hours of portable operation and a bright digital readout. Diagnostic quality waveforms can be captured and reviewed through the Synergy®-E Software. In addition, an easy-to-understand interface incorporates universally recognized symbols that help reduce potential language barriers.

SmartMonitor 2 PS, however, is not just for infant use. Additional pediatric and adult clearances make the SmartMonitor 2 PS ideal for documenting patient response to conscious sedation, patient-controlled analgesia and general floor monitoring.

The SmartMonitor 2 PS offers:

- Superior waveform quality
- 4 MB of memory capacity
- 15 hours of portable operation
- Universal symbol interface
- Lightweight design
- Hospital or homecare selectable settings
- Backup battery
- Dual alarm
- Durable metal patient connector



Synergy-E Software provides graphic displays for reviewing event data.



SPECIFICATIONS

Dimensions	2.25 in x 7.25 in x 9 in (5.7 cm x 18.4 cm x 22.9 cm)
Weight	3.0 lbs (1.35 kg)
Internal Channels	QRS (ECG), Respiration, SpO ₂ , Pulse, Plethysmograph
External Channels	Four Auxiliary Channels

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	QUANTITY
1014557	SmartMonitor 2 Professional Series	1/each
1030271	SmartMonitor 2 Professional Series with Hospital Alarm	1/each
1030193	Replacement Metal Patient Cable (Hospital 0.060")	1/each
1030192	Replacement Metal Patient Cable (Homecare 0.080")	1/each
1014558	Replacement Plastic Patient Cable (Hospital 0.060")	1/each
H4005	Replacement Plastic Patient Cable (Homecare 0.080")	1/each
H5870	Masimo Oximeter Patient Cable	1/each
1025721	Masimo LNOP Neo-L Neonatal Disposable Sensor	20/box
1025723	Masimo LNOP Infant-L Infant Disposable Sensor	20/box
1025725	Masimo LNOP NeoPt-L Preemie Disposable Sensor	20/box
1015662	SmartTrace™ ECG Neonatal/Pediatric Electrodes (Homecare 0.080") (2 electrodes per set; 20 sets per case)	20/case
1015663	SmartTrace™ ECG Neonatal/Pediatric Electrodes (Hospital 0.060") (3 electrodes per set; 10 sets per case)	10/case
1022895	USB to Serial Adapter	1/each
1022029	Synergy®-E Software V3.10	1/each
1028970	SmartMonitor 2 Professional Series Light with Standard Alarm (85dB)	1/each
1038140	SmartMonitor 2 Professional Series Light with Hospital Alarm (70dB)	1/each



Customer Service: 800-345-6443 or 724-387-4000
 Respirationics Europe, Africa, Middle East: +33-1-47-52-30-00
 Respirationics Asia Pacific: +852-3194-2280
www.childmed.com

Respirationics, Synergy and SmartMonitor are registered trademarks and the Children's Medical Ventures heart design logo and SmartTrace are trademarks of Respirationics, Inc. its affiliates.

Masimo SET is a registered trademark of Masimo Corporation.

©2008 Respirationics, Inc. and its affiliates. All rights reserved.